
Development of low background CsI(Tl) and NaI(Tl) crystals for WIMP search

Hyun Su Lee

Ewha Womans University

On behalf of the KIMS Collaboration

Low Radioactivity Techniques 2015

KIMS (Korea Invisible Mass Search)

- Dark matter search at **Yangyang underground laboratory**
 - ❖ May have new laboratory in Samcheok
- Funded by National Research Foundation of Korea (2000)
 - ❖ Dark matter (DM) search with CsI(Tl) crystals (KIMS-CsI)
- Establishing the Center for Underground Physics (CUP) in the Institute of Basic Science (IBS) (2013)
 - ❖ Upgrade of KIMS-CsI
 - ❖ DM searches with NaI(Tl) crystals (KIMS-NaI)
 - ❖ DM searches with low temperature detector (KIMS-LT)

Yangyang Underground Laboratory

(Upper Dam)

Korea Middleland Power Co.
Yangyang Pumped Storage Power Plant

700 m



(Power Plant)



Minimum depth : 700 m
Access to the lab by car (~2km)

Internal background of CsI(Tl) crystals at beginning


The lowest background CsI(Tl) crystal in the market




Nucl. Instrum. Meth. A 500 (2003) 337

^{137}Cs reduction

- Water is main source




Nucl. Instrum. Meth. A 552 (2005) 456




- With **purified water**, we can reduce the internal background

Internal background of CsI(Tl) crystal with “pure water”



Nucl. Instrum. Meth. A 571 (2007) 644

- One test and Two full sized crystals were developed
- Physics results with pilot run




Phys. Rev. Lett. 99 (2007) 091301


Phys. Lett. B 633 (2006) 201

CsI(Tl) crystals with “ultra pure” water

- Using “ultra-pure” water
- Growing **twelve full sized crystals**




Multiple –hit backgrounds




- **Three years stable operation**
(2009-2012)

Physics results with a twelve CsI crystal array

SI cross section limit



Annual modulation analysis




S.C. Kim et al., PRL 108 181301 (2012)

H.S. Lee et al., PRD 90 052006 (2014)


Used 1/3 data set

Full data set result soon will be out

< SI WIMP-nucleon Cross Section >



KIMS-NaI experiment



- To confirm DAMA annual modulation signature
 - ❖ CsI is not enough for WIMP-Na interaction
 - ❖ Same NaI crystal for the same annual modulation signature
- Need to develop ultra-pure NaI(Tl) crystals
 - ❖ Goal is less than DAMA background ($\sim 1 \text{ dru} = 1 \text{ counts/keV/kg/day}$)
 - ❖ 200 kg× 3 years data will prove DAMA signature without any ambiguity

KIMS-NaI crystals

- Development of low background NaI(Tl) crystals

	Nal-001	Nal-002	Nal-003	Nal-004	Nal-005	Nal-006
Mass	8.26 kg	9.15 kg	3.35kg	3.35kg	9.16 kg	11.44 kg
Powder	AS	AS	SA-AG	SA-CG	AS	SA-CG
Crystal	AS	AS	AS	AS	AS	BH
Arrive	2013.9	2014.1	2014.8	2014.8	2014.11	2014.12

K.W.Kim et al., Astropart. Phys. 62, 249 (2015)

Glossary


AS = Alpha Spectra Inc (US company)

SA-AG = Sigma Adrich, Astro-grade (less K40)



SA-CG = Sigma Adrich, Crystal-grade

BH = Beijing Hamamatsu (China)

KIMS-NaI detector module



- Hamamatsu R12669 PMTs are attached
 - ❖ Supposed same PMTs with recently upgraded DAMA PMTs
 - ❖ ~35% quantum efficiency at 420nm
- Light Yield: ~15 photoelectrons/keV
 - ❖ Consistent with ANAIS-25
- Data taking
 - ❖ 400MHz Flash ADC (Notice Korea)
 - ❑ Flexible trigger logic with FPGA
 - ❑ Trigger condition: 1 PE/PMT within 100ns



**~200 ns
decay time**


PMT background

- We can well identify PMT noise events using same parameter with DAMA

NIMA 592, 297 (2008)


DAMA

$2 \text{ keV} \leq E \leq 4 \text{ keV}$




$$X_1 = \frac{\text{Area (from 100 to 600 ns)}}{\text{Area (from 0 to 600 ns)}}$$

$$X_2 = \frac{\text{Area (from 0 to 50 ns)}}{\text{Area (from 0 to 600 ns)}}$$




Astropart. Phys. 62, 249 (2015)

KIMS-NaI




Remained noise?


Single Hit events



Multiple Hit events




$$Asymmetry = \frac{(Q_{PMT0} - Q_{PMT1})}{(Q_{PMT0} + Q_{PMT1})}$$





We should apply additional cuts to reject such noisy events

Intrinsic Background – ^{40}K



DAMA crystals :
10~20 ppb

Coincidence signals



Geant4 simulation
on efficiency

	NaI-001	NaI-002	NaI-003	NaI-004	NaI-005	NaI-006
K (ppb)	41.4 ± 3.0	49.3 ± 2.4	25.3 ± 2.4	> 110	40.1 ± 4.2	> 150

Intrinsic Background – ^{40}K


	Nal-001	Nal-002	Nal-003	Nal-004	Nal-005	Nal-006
K (ppb)	41.4 ± 3.0	49.3 ± 2.4	25.3 ± 2.4	> 110	40.1 ± 4.2	> 150
Powder	AS	AS	SA-AG	SA-CG	AS	SA-CG
K (powder)	?	?	25.07	~200	?	~200

- **Nal-003** used Sigma Adrich **Astro-Grade** (SA-AG) powder **25 ppb**
- Nal-004 & Nal-006 used **Crystal-Grade** (SA-CG) powder
- All Alpha Spectra prepared powders (Nal-001, Nal-002, Nal-005)
~ 40 ppb levels

Nal powder is a key of K contamination

- Will grow one another crystal using different batch of SA-AG powder
❖ ~10 ppb from ICP-MS by Sigma Aldrich
- **K measurement technique** from powder is under development
- R&D to **reduce K in Nal powder** was just started!!


Intrinsic Background (alpha analysis)




Radionuclei	NaI-001 [mBq/kg]	NaI-002 [mBq/kg]
^{238}U (^{214}Bi)	<0.007	<0.001
^{228}Th (^{216}Po)	<0.012	0.002 ± 0.001
^{210}Po	3.28 ± 0.02	1.76 ± 0.01
Total alphas	3.29 ± 0.02	1.77 ± 0.01

^{238}U and ^{228}Th contaminations were very small!!


^{232}Th chain



^{238}U chain



Intrinsic Background – ^{210}Pb




- Most of alphas are coming from ^{210}Po
- Broken equilibrium (Rn-222 contamination) during powder or crystallization

	Nal-001	Nal-002	Nal-003	Nal-004	Nal-005	Nal-006
Total alpha (mBq/kg)	3.29 ± 0.01	1.77 ± 0.01	2.29 ± 0.03	-	0.47 ± 0.01	1.81 ± 0.02
Powder	AS	AS	SA-AG	SA-CG	AS	SA-CG


- **Nal-003, Nal-004, and Nal-005** had better treatment on air (Rn) contact during crystal growing
- **Nal-005** had better treatment of chemical process on powder

Intrinsic Background – ^{210}Pb


Nal-002




Nal-003



Nal-005




Nal-006



Intrinsic Background – ^{210}Pb

- Main Source.. Powder production or crystal growing?
 - ❖ Should know ^{210}Pb in powder
 - ❖ Alpha rate measurements at powder
 - ❑ Buy High sensitive alpha counters
 - ❑ Developing ZnS detector
 - ❖ 46 keV gamma measurement
 - ❑ Considering Well-type Ge or Si detectors



Kamland-pico demonstrates a reduction of ^{210}Pb in NaI powder with ion-exchange resin (arXiv:1407.3542): $50 \mu\text{Bq/kg}$ in crystal

We just started ^{210}Pb reduction in powder with various resins
Co-worked with Russian chemists

Nal Powder purification

Sorbent	Matrix	Application
A	Inorganic, manganese dioxide	K ⁺ , Pb ²⁺
B	Inorganic, zinc ferrocyanide	K ⁺
C	Organic macroporous anion exchange copolymer	Pb ²⁺
D	Organic macroporous chelating copolymer	Pb ²⁺
E	Organic, Resorcinol-formaldehyde	K ⁺

R&D for K and Pb purification

- ❖ First : Simple shaking for 24hrs with resins
 - ❑ Just done a week ago
- ❖ Next : Column chromatography

Shaking NaI solution with various resins





Preliminary result of purification

- Measurement with ICP-MS for K and Pb
 - ❖ We made NaI solution with some amounts of K and Pb

Resins	³⁹ K (ppb)			²⁰⁸ Pb (ppb)			Applications	
	Bef.	Aft.	DF	Bef.	Aft.	DF		
A	223	254	0.9	1275	19	67	K ⁺ , Pb ²⁺	
B	463	375	1.2	1693	1589	1.1	K ⁺	
C	404	313	1.3	1653	5	330	Pb ²⁺	
D	352	312	1.1	1619	33	49	Pb ²⁺	
E	378	340	1.1	1661	898	1.8	K ⁺	

- Reduction of **Pb** seems to **work very well**
- However, **K** seems to be **difficult**
- We will have much more experiences

Internal background




Simulation for external and cosmogenic activation is under development

- We achieved ~ 2 dru at 6keV
- ~ 0.5 dru caused by internal backgrounds (^{210}Pb)
 - ❖ Additional ~ 0.7 dru at 2-4 keV due to ^{40}K
- ~ 1.5 dru caused by **external** (+cosmogenic)
- **Goal : less than 1 dru (both less than 0.5 dru)**

External background reduction

- We prepared liquid scintillator active veto system
 - ❖ **Veto efficiencies** for sources from PMT radioisotopes (U, Th, K) were greater than **80%** at low energy (0-10 keV)

Liquid Scintillator



NaI crystal

**Less than 0.5 dru
from external might
be possible**

	Feb/15	Mar/15	Apr/15	May/15
Shielding Box	Design & Construction			
LS base on LAB	Propagation & Insert			
Run			Run&Analysis	

Prospect of NaI(Tl) crystal development



- ^{40}K less than 10ppb
 - ❖ Astro-grade powders from Sigma-Adrich with ^{40}K less than 10ppb or our own purification (**<0.2 dru** @2-4 keV)
 - ❖ ~factor two reduction with liquid scintillator veto system
- ^{210}Pb less than 0.2 mBq/kg
 - ❖ Already 0.5 mBq/kg
 - ❖ At least factor 2 reduction with powder purification (**<0.25 dru** @2keV)
- External
 - ❖ **less than 0.5 dru** with liquid scintillator veto system
- **Goal : 200kg NaI detectors (less than 1 dru background) within 2 years**

Conclusion

- KIMS collaboration successfully grew the low background CsI(Tl) crystals for the WIMP search
 - ❖ Have been contributed for understanding of DM nature during late decade
- **KIMS-NaI** collaboration start to grow the ultra-low background NaI(Tl) crystals for the WIMP search
 - ❖ Current ~2dru @ 6keV
 - ❖ Goal ~1dru @ 2keV (might be possible within one year)
 - ❖ We will **confirm the DAMA** signature within a couple of years


Backup

Underground laboratory



- **Underground space for the KIMS-Nal was already prepared in Y2L**
 - ❖ 10,000 class clean room
 - ❖ Rn free air will be supplied
 - ❖ Size of detector room ~50 m²

Preliminary design of KIMS-NaI



- 200 kg NaI
 - ❖ 30 cm liquid scintillator veto (Active veto)
 - ❖ 6 cm Cu
 - ❖ 20 cm lead
 - ❖ Muon veto (3 cm plastic scintillator)
 - ❖ 30 cm Polyethylen